FICOLL PROTOCOL FOR WBC EXTRACTION


1 Obtaining MONONUCLEAR cells of peripheral blood with LEUCOSEP tube

Optional: put together the content of all them tubes of extraction of blood in a single tube falcon of 50 ml to homogenize the blood.


1. Add 3,2 ml of Ficoll to each labeled tube Leucosep 10 ml (empty tubes). It depends of the volume of LEUSEP tube. The ficoll must reach the filter.

2. Centrifuge the Leucosep tube containing the ficoll to 1000xg for 1 min. Remove the excess Ficoll with pipette or decantation.
3. Centrifuge tubes of blood at 1500 g for 15 min
4. Vacuum plasma without catch mononuclear cells that form a white layer above the erythrocytes and granulocytes and discard.
5. Sucking 2 ml of the mononuclear cell layer and put it in a tube falcon 15 ml containing 4 ml of PBS. Invert the tube to mix the solution well.
6. Add the layer of mononuclear cells with PBS or blood to the Leucosep tube of point 1.

7. Centrifuge tubes to 800 g for 15-30 min. to 18-25 ° C without brake


From here continue at point 3.1


2 Obtaining MONONUCLEAR cells of peripheral blood without LEUCOSEP tube

Optional: Gather the contents of all the tubes of blood in a single 50ml falcon tube to mix the blood.


1. Add a volume of Ficoll for each two volumes of blood (Ficoll: blood 1:2)
2. Centrifuge tubes of blood to 1500g for 15 min.
3. Vacuum plasma without catch mononuclear cells that form a white layer above the erythrocytes and granulocytes and discard. Add PBS in amount equal to the withdrawn plasma volume.
4. Mix the blood with the PBS and pass to the next step.
5. Sucking 2 ml of the mononuclear cell layer, pour it into a Falcon 15 ml tube containing 4 ml of PBS. Invert the tube to mix the solution well.
6. Add the layer of mononuclear cells with PBS or as horizontal as possible (drop to drop in the wall of the tube, slowly), the blood on the ficoll with tube so the ficoll and blood do not mix.

7. Centrifuge the tube at 500 g for 10 min at 18-25 ° C without brake. (type of rotor: swing-arm rotor)

Web to calculate your RPM/G-force (<http://www.endmemo.com/bio/grpm.php>)

From here continue at point 3.1

3. Obtaining of MONONUCLEAR cells from peripheral blood

1. Remove the White halo of CMSPs carefully with a Pasteur pipette and transfer to a sterile tube. Make up the volume with PBS.
2. Centrifuge the tube at 500 g for 10 min at 18-25 ° c
3. Remove tubes supernatant by decantation taking care not to break the cell button. Resuspend the pellet in PBS.
4. Extract a sample of 10-25 µl and carry out a count.
5. Centrifuge the tube at 500 g for 10 min at 18-25 ° c
6. Discard the supernatant and resuspend the pellet in half-full:
 - For cultivation (RPMI + 10% FBS)
 - To freeze (SBF+7, 5 - 10% DMSO)
7. Perform the appropriate dilution to get a desired concentration:
 - For freezing 7 - 15 x 10⁶ cel. ml.
8. freeze or put into cultivation.

FICOLL PROTOCOL FOR WBC EXTRACTION(Spanish)

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